



11/2018

AGATA Data Analysis User's guide

Local Level Processing

AGATA Data Analysis Team

This document provides a guide to help the users analyzing the AGATA data produced at the local level processing i.e. before any building of events. It includes traces / energy / time calibrations, alignments, cross talk corrections and any other corrections to improve the quality of the data. Criteria for bad events rejection are also highlighted.

The last version of this document, provided by the AGATA Data Analysis Team, can be found on ATRIUM (<https://atrium.in2p3.fr/>). Just click on WORSPACE (ARBORESCENCE) and follow the tree on the left side (see the following picture) to reach the DataProcessing entry point.

Document, utilisateur, groupe

HOME ARBORESCENCE RECHERCHE TICKET SUPPORT ME CONNECTER Guest

Atrium Zone de Diffusion Grand Public Projet AGATA DataProcessing

DataProcessing

Documents related to online/offline processing of the Data

Contenu Résumé Historique

Documents / page: 30

Type	Titre	Atrium ID	Créé le	Auteur	Modifié le	Modifié par	Vrs. Obs.	Etat
Document	WatcherUsersGuide	ATRIUM-3982	5 mars 2015	Olivier STEZOWSKI	21 mai 2015	stezow@ipnl.in2p3.fr	3.0	Validé
Document	GwUserGuide.odt	ATRIUM-3884	27 févr. 2015	Olivier STEZOWSKI	27 févr. 2015	stezow@ipnl.in2p3.fr	1.0	Validé
Document	Cookbook.odt	ATRIUM-3883	27 févr. 2015	Olivier STEZOWSKI	27 févr. 2015	stezow@ipnl.in2p3.fr	1.0	Validé

Contents

Some general statements.....	3
AGATA Forum.....	3
NARVAL, actors and the AGAPRO package.....	3
The AGATA local level Processing.....	3
Replay of Data / Directories organization.....	5
Replay of Data / Emulators.....	6
Narval (and Narval standalone).....	7
Femul.....	7
The different actors processing the data flow.....	11
The Producer actor.....	11
The Preprocessing actor.....	12
Energy calibration.....	12
Pole-zero / shaping-time adjustments.....	21
Crosstalk corrections.....	22
Time alignment.....	29
The Pulse Shape Analysis (PSA) actor.....	32
The PostPSA actor.....	32
neutron damage corrections.....	32
Final energy re-calibrations (with offsets).....	36
Force Segments to Core.....	39
Global time alignments.....	39

People involved in this document:

D. Bazzacco, A. Boston, E. Clement, N. Dosmes, J. Dudouet, A. Gadea, F. Holloway, L. Hongjie, A. Korichi, N. Lalovic, E. Legay, J. Ljungvall, C. Michelagnoli, R. Perez, D. Ralet, M. Siciliano and O. Stezowski

NOTE: Please for any comment/question/suggestion contact agata@ipnl.in2p3.fr

Some general statements

AGATA Forum

The AGATA Forum <http://agata.in2p3.fr/forum/> is a dedicated place for any question, issue, or discussion concerning the AGATA data analysis.


NARVAL, actors and the AGAPRO package

The AGATA DAQ box is based on NARVAL [Ref, see also the section '**the ADF component**' in **the GammaWare User's Guide**]. As a consequence of such a framework, the AGATA data is processed through consecutive calls of ACTORS, each actor being in charge to delivering useful data for the next one in the chain. Even if NARVAL is itself written in ADA, it can bind other languages (C, C++, ...). However, the choice of the AGATA collaboration has been to develop algorithms in the C++ language. The most relevant actors are part of a package names AGAPRO [see **How to install the AGAPRO package in the Cookbook**].


The AGATA local level Processing

At the local level processing, three actors, also called filters in NARVAL's terminology, process the data flow. Those actors, Preprocessing, PSA (Pulse Shape Analysis) and PostPSA¹ are in charge respectively to prepare traces (energy calibration / time alignments), to determine hits from traces and to apply additional corrections/filters to prepare the data to be merged. At this stage one achieves the global level processing which is mainly based on the gamma-ray tracking. Such a chain is illustrated in the following picture:

¹ Depending on the data flow processed, this filter may not be mandatory.



Before any processing, an actor is initialized first using a configuration file which has in principle the *.conf* extension. Such a file contains all the parameters required by the algorithm to work properly. This is illustrated in the following picture:



NOTE:

To manage all the configuration files needed by all the actors, a python script, called *gen_conf.py*, is provided.

In addition, to determine the best parameters for a given data set, specific works (or actions) are required which are based on spectra and/or on more complex inputs (traces, events, energy) at different levels in the processing of the data flow. Several mechanisms are available to do this. Traces/events can be read from the data stream (in the ADF format) out of any actors and thus can be used to build spectra. This is the way the Watchers work. For historical reasons, spectra, traces/events can be produced also directly by the actor itself² and saved on disk together with the ADF files.

The spectra produced by the actors are saved in a specific format that can be read by TkT or by the GammaWare package. Specific treatments are then respectively realized through dedicated command line programs, available in the AGAPRO package, or through graphical user's interfaces in the GW/Watchers.

Replay of Data / Directories organization

The data produced by an experiment are stored in a single (unix) directory itself composed of several sub-directories one per run. For each run, one can find the data produced (*Data* sub-

² Depending on the configuration of the processing chain at compilation and from the different .conf files

directory) at different level in the data flow as well as the files used to configure (*Conf* sub-directory) the different actors.

/agatadisks/eXXX/eXXX/run_0001.dat.14-04-17_11h12m53s

- |-- **Conf** → Contains the configuration files of each actor
- |-- **Data** → Contains the raw data (*.adf and *.cdat)
- `-- **gen_conf.py** → Python script used for actors settings (more details in the following)

The two sub-directories are themselves divided in sub-directories, one per Germanium crystal for the local level processing and several ones for the ancillaries and the global level processing.

/agatadisks/eXXX/eXXX/run_0001.dat.14-04-17_11h12m53s

- |-- **Conf**
 - |-- **00A**
 - ...
 - |-- **14C**
 - |-- **Builder**
 - |-- **Global**
 - |-- **Merger**
- `-- **Data**
 - |-- **00A**
 - ...
 - |-- **14C**
 - |-- **Builder**
 - |-- **Global**
 - |-- **Merger**

Here is also a snapshot of a typical content of the bottom directories.

- |-- **14C**
 - |-- **BasicAFC.conf**
 - |-- **BasicAFP.conf**
 - |-- **CrystalProducer.conf**
 - |-- **CrystalProducerATCA.conf**
 - |-- **CrystalProducerGGP.conf**
 - |-- **PSAFilter.conf**
 - |-- **PostPSAFilter.conf**
 - |-- **PreprocessingFilter.conf**
 - |-- **PreprocessingFilterPSA.conf**
- |-- **Builder**
 - |-- **BasicAFC.conf**
 - |-- **CrystalPositionLookUpTable**
 - |-- **EventBuilder.conf**
 - `-- **TrackingFilter.conf**
- |-- **Merger**
 - |-- **BasicAFC.conf**
 - |-- **CrystalPositionLookUpTable**
 - |-- **EventMerger.conf**
 - `-- **TrackingFilter.conf**

Replay of Data / Emulators

In the process to optimize the parameters to get the best quality as possible for your data, it is required to apply one or several times the same algorithm, or actor, using different configuration

files. This could be done using NARVAL itself, on site or out site. Another possibility is to use what is called emulators i.e. other frameworks in charge of organizing / running the processing of data by many different actors.

Narval (and Narval standalone)

How to replay with NARVAL on site / out site ... To be written

Femul

FEMUL (Flat EMULator) is able to run complex topologies including DISPATCHER i.e. actors having several input lines and one output line. There are two such dispatchers, *EventBuilder* and *EventMerger*: the former being in charge of building AGATA events (crystals in coincidence from the PSA_XXXX.adf and resulting in Builder_XXXX.adf files) and the later allows one to build coincidences with ancillaries (resulting in merged_XXXX.adf files). All actors are part of the AGAPRO package. The procedure to get/compile/install this program is given in the Cookbook [see **How to install the femul emulator in the Cookbook**].

The topology (list of detectors and actors to apply on the data flow) is built using an ascii “Topology.conf” file and the way to run is:

```
femul Topology.conf3
```

Depending on the task you would like to perform, you may need to adapt the topology given in *Topology.conf*. The list of the most commonly used actors is:

Producer actors (to start a NARVAL chain):

CrystalProducerATCA → To start a replay from traces (*.cdat) files

BasicAFP → To start a replay from ADF files (Basic Agata Frame Producer)

Filter actors (to apply calibrations/algorithms on the data flow):

PreprocessingFilterPSA → To apply energy and time calibrations/corrections

PSAFilter → To execute the Pulse Shape Analysis algorithm

PostPSAFilter → For final corrections after PSA (neutron damage, re-calibrations...)

TrackingFilterOFT → To process the tracking algorithm

Dispatcher actors (to merge several NARVAL chains in one):

EventBuilder → To build a AGATA event within a specific timing window

EventMerger → To build a global event between AGATA and an ancillary detector in a specific timing window

Consumer actors (to end a NARVAL chain):

BasicAFC → To end a replay and write ADF files on disk (Basic Agata Frame Consumer)

³ femul -h to get some help from the command line

None → To end a replay without writing the output files (used when only spectra files are used in calibration procedures)

Here are three different examples of such topology file:

The first one shows how to read traces (fevent_mezzdata.cdat.*) from many crystals, run up to the PSA and then dump hits in ADF files. (eg: PSA_XXXX.adf)

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 4    CRY
Producer   CrystalProducerATCA
Filter     PreprocessingFilterPSA
Filter     PSAFilter
Consumer   BasicAFC
ENDLOOP
```

The second example shows how to read psa hits (psa_*.adf) from many crystals, perform some post PSA task, build AGATA events, apply tracking and then dump tracked gamma-rays in ADF files (Tracked_*.adf)

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 3    CRY
Producer   BasicAFP
Filter     PostPSAFilter
Dispatcher EventBuilder
ENDLOOP

Chain 3    Global/
Builder    EventBuilder
Filter     TrackingFilter
Consumer   BasicAFC
```

The last example shows how to read psa hits (psa_*.adf) from many crystals, perform post PSA, build AGATA events, merge them with some ancillary data before applying tracking and then dump tracked gamma-rays in an ADF file.

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 3    CRY
Producer   BasicAFP
Filter     PostPSAFilter
Dispatcher EventBuilder
ENDLOOP

Chain 3    Builder/
Builder    EventBuilder
Consumer   BasicAFC
Dispatcher EventMerger

Chain 2    ancillary/
Producer   BasicAFP
Dispatcher EventMerger

Chain 3    Merger/
Builder    EventMerger
```


Filter	TrackingFilterOFT
Consumer	BasicAFC

To replay data, one should create a top destination directory containing the same structure as the one produced online. The *Conf* sub-directory and the *gen_conf.py* should be copied from the initial one, and the original Data should be linked to avoid copies of large amount of data. The *gen_conf.py* is then used to adapt the configuration files to your current directory. It will also produce an *Out* sub-directory where the files produced by the replay will be stored. Here is an example to prepare a replay folder for run_xxxx, where the raw data are located in /agatadisks/:

```
$ mkdir Replay
$ cd Replay
$ mkdir run_xxxx
$ cd run_xxxx
$ cp -r /agatadisks/run_xxxx.../Conf ./
$ cp /agatadisks/run_xxxx.../gen_conf.py ./
(Here edit the gen_conf.py as explain in the following)
$ ln -s /agatadisks/run_xxxx.../Data Data
$ python gen_conf.py
```

The command: `python gen_conf.py -h` show a help message

When copying the *gen_conf.py*, take care of changing **ONLINE** into **OFFLINE** and **NARVAL** into **femul**.

Here is the beginning of the *gen_conf.py* script showing some parameters to be changed for the replay you would like to perform.

```
Script to generate the configuration files for the replay of AGATA data using the
multi-process distributed system Narval or the single-process emulator femul.
The script is divided in a few sections, some of which are specific to the actual
analysis (and therefore are likely to be changed by the user) while some others are
normally not to be touched:
0) Type of analysis and replacement macros (in the style of the shell) used to parametrize
the commands listed in 2).
1) The structure of the actual analysis is defined by the variables PROGTYPE and CONFTYPE and
by the dictionaries Topology and Actors. The dictionary ExtraFiles contains a list of files,
which are needed by the analysis but are not generated by this script (e.g. calibrations,
mappings, ...); these files can be copied from a previous analysis if the script is started
with the option -o or --old (python gen_conf.py -h to get the list of accepted options)
2) The command lines to be written in the conf files of the actors defined in Actors{}.
Uncomment/comment/modify the command lines and their parameters
3) A small database defining the position and the PSA signal basis of the germanium crystals.
This part should not need to be changed.

To get a list of command line arguments, launch the script as:  gen_conf.py -h
To get a list of keywords accepted by the vaious actors:      femul -k
"""
#####
##### 0 Type of analysis and replacement symbols #####
#####
```

```

PROGTYPE='femul'      # NARVAL or femul      (to choose between os.getcwd() and '' for CWD)
CONFTYPE='OFFLINE'    # ONLINE or OFFLINE (used just to exclude the ReadDataDir line in the
Producers)

MACROS={              # various replacements for symbols defined in 2).
'$CONFDIR'           : 'Conf',              # this will be prefixed by CWD/
'$READDIR'           : 'Data',              # this will be prefixed by CWD/; if ONLINE this will
not be written
'$SAVEDIR'           : 'Out',               # this will be prefixed by CWD/; if ONLINE this will
be replaced by $READDIR
'$ANCILLARY'         : 'Ancillary',          # this will be prefixed by CWD/
'$GLOBAL'            : 'Global',            # this will be prefixed by CWD/
'$PSABASE'           : '/agatadisks/bases/ADL', # standard place at LNL/Linux
'$CRYSTAL_ID'        : "",                  # the actual value is defined in GeDataBase
'$SIGNAL_BASIS'      : "",                  # the actual value is defined in GeDataBase
'$CRYSTAL'           : "",                  # the actual value taken from Topology['CRYSTAL']
}

#####
##### 1 Structure of analysis #####
#####

Topology={           # The directories to be generated in Conf, Data and Out
'CRYSTAL'           : "00A 00B 00C 01A 01B 01C 02A 02B 02C 03A 03B 03C 04A 04B 04C 05A 05B 05C",
'ANCILLARY'         : "Ancillary",
'BUILDER'           : "Builder",
'MERGER'            : "Merger",
}

# The name of the used actors must correspond to one of the tuples defined in the following
section.
# This requirement creates a problem for BasicAFP and BasicAFC when they are used in chains of
different type
# (e.g. after PSA and after Tracking) and one wants to define chain-specific names for their
input/output files.
# The solution is to suffix the name of the chain-type (e.g. _CRYSTAL or _GLOBAL or any other), to
the defining tuple.
# This suffix will be silently removed from the actual name of the generated configuration files.

Actors={             # These are the xxxx.conf files to be generated
'CRYSTAL'           : "CrystalProducer BasicAFP_CRYSTAL PreprocessingFilter PSAFilter PostPSAFilter
BasicAFC_CRYSTAL",
'ANCILLARY'         : "BasicAFP_ANCYLLARY AncillaryFilter BasicAFC_ANCYLLARY",
'BUILDER'           : "BasicAFP_BUILDER EventBuilder BasicAFC_BUILDER",
'MERGER'            : "EventMerger TrackingFilter BasicAFC_MERGER",
}

ExtraFiles={         # If not already present, these files can be copied from a directory specified in the
command line
'CRYSTAL'           : "CrystalProducerATCA.conf PreprocessingFilterPSA.conf xinv_1325-1340.cal xdir_1325-
1340.cal",
'ANCILLARY'         : "AncillaryFilter.conf",
'BUILDER'           : "CrystalPositionLookUpTable",
'MERGER'            : "CrystalPositionLookUpTable",
}

#####
##### 2 Tuples specifying the content of the configuration files #####
#####

```

In case of replay of the PSA actor, the ADL bases need to be downloaded from the AGATA DAQ box or from the Cologne web-site : <https://www.ikp.uni-koeln.de/research/agata/index.php?show=download> (contact agata{at}ipnl.in2p3.fr for download rights). The path to the ADL bases needs to be updated in the *gen_conf.py* script (in the `$PSABASE` variable).

The CrystalProducer actor needs also to be modified to specify to use the traces files as input. For this, uncomment the lines “InputDataFile” and “AllInputFiles”:

```
#####  
CrystalProducer=(  
  "ActualClass          CrystalProducerATCA",  
  "CrystalID            $CRYSTAL_ID",  
  "ReadDataDir          $READDIR/$CRYSTAL",  
  "SaveDataDir          $SAVEDIR/$CRYSTAL",  
  "TraceLength          100",  
  "WriteDataMask         0",  
  "InputDataFile         event_mezzdata.cdat",  
  "AllInputFiles",  
)
```

Finally, make sure that the option “NoMultiHist” is commented in the different actors of the gen_conf.py during the calibration procedure. This allow to produce the spectra files used in the different steps of this document.

The different actors processing the data flow

The Producer actor

This actor is in charge to start the NARVAL chain. In online mode, it gets the traces from the front end electronics and sends it to the data flow. In offline mode, it reads input files (adf or traces files)

and sends them to the NARVAL chain.

The Preprocessing actor

This actor is in charge of preparing the data for the PSA algorithm. In order to do that job, two files which are named **PreprocessingFilter.conf** and **PreprocessingFilterPSA.conf** should contain the required information.

As any detector, any AGATA capsule should be calibrated in energy. Since it is also highly electrically segmented, cross talks have a significant importance and must be corrected. For details, see for instance these publications for such effects: B.Bruyneel et al., NIMA 608 (2009)

Be aware, there might be time delays between the different segments and the core signal. Any PSA code requires having all the signals from one capsule perfectly aligned and thus time alignment are also performed.

Here are snapshots of the two required configuration files.

For the first one (“**PreprocessingFilter.conf**”):

ActualClass	PreprocessingFilterPSA
SaveDataDir	/global/path/to/Replay/Folder/run_xxxx/Out/00A
EnergyGain	4
XtalkFile	xinv_1325-1340.cal
WriteTraces	100

And the second one (“**PreprocessingFilterPSA.conf**”):

segm	0	4600	500	0.300652	15	17.220
segm	1	4600	500	0.313630	15	17.873
...						
core	0	4390	500	0.482203	0	20.
core	1	4390	500	1.735468	0	20.
tntf	2097152					

Note: the content may be slightly different depending on the version of the PreprocessingFilter used. For the second file, here is the meaning of the different columns:

type (segm/core) segment/core ID TFall TRise egain emink tmove

TFall Pole zero / decay time [timestamp units]

Trise Shaping time / Risetime [timestamp units]

egain slope of energy calibration (no offset at the preprocessing level) [keV/channel]

emink energy threshold [keV]

tmove Shift to align in time the different segments [ns]

Energy calibration

Explanations and goals

In digitizers, signals are processed to extract the amplitude using a trapezoidal filter. This value is written from the beginning of the chain into the data flow. This is the value which is used to calibrate the detectors in energy. Concerning the 36 segments, because of the various cross talks, the calibration in energy is done using events in which only one and only one segment in a given crystal

has fired.

The obtained calibration coefficients are to be set in the 5th col of the *PreprocessingFilterPSA.conf* file.

Note: All the different steps of the following procedures needs to be done for each crystal. The environment variable `$CryId` will be used in the following to refer to the current detector (ex: `CryId=00A`). We suppose for the given shell command that you are working in the Conf folder of the current detector.

Tools available

**** command line programs ****

To use those tools, spectra directly produced by the actors at running time are required. In particular the ones contained in this library⁴:

`Data/$CryId/Prod__4-38-32768-UI__Ampli.spec`

From the file name, we can deduce that this spectra file contains 4 libraries of 38 spectra written in 32768 unsigned integer bins, containing the amplitude spectra of segments and cores. The first library contains the raw amplitude from the trapezoidal filter. The second library contains events of segment multiplicity=1, with amplitude larger than the threshold and scaled by the EnergyGain (see conf file) factor. The calibration needs to be done on spectra from the second library. The two other libraries are not used in the energy calibration but can be necessary for x-talk corrections.

From these spectra files, a C program, called *RecalEnergy*⁵, is able to find peaks and thus calibrate the different channels (core and segments) for one crystal. Here is the way to use it, with some options⁶:

```
$ cd Conf/$CryId # If not already done, go in the working directory
$ RecalEnergy -spe ../../Data/$CryId/Prod__4-38-32768-UI__Ampli.spec -sub 38 -num 38
-gain 2 -60Co |tee recal.out
```

Options:

- `-spe filename` : name of spectrum to analyze (mandatory)
- `-sub nn` : analysis starts from spectrum nn (`-sub 38` to analyze the second library)
- `-num nn` : number of spectra to analyze (`-num 38` for 36 segments + core high and low gain)
- `-lim` : limit the search to this range in channels (use if necessary)
- `-dwa` : default fwhm and minimum amplitude for the peaksearch
- `-60Co` : define the source (default is 60Co, more than one source is allowed)
- `-gain val` : scaling factor for the slope

The gain factor needs to be coherent with the gain applied to the spectra, this information can be found in the file *CrystalProducer.conf* : `ProjeM1 10 2` (threshold and gain)

A way to apply this for several crystals `$DetList` is the following:

Note: You need to be in the global folder (where Out, Data, Conf are located)

⁴ These spectra are produced by the CrystalProducer actor.

⁵ See zPrograms directory in the agapro package

⁶ For a full list of options, `RecalEnergy -h`

```
$ for i in $DetList ; do RecalEnergy -spe Data/$i/Prod__4-38-32768-UI__Ampli.spec -sub
38 -num 38 -gain 2 -60Co |tee Conf/$i/recal.out ; done
```

In the following, the commands will be written for only one detector, but this loop can be use to obtain automatic procedures on a list of detectors.

Check carefully (rEnergy, FM05, Chi2...) the output of each file to be sure that the fits are good. If not, reprocess the concerned files using more restrictive parameters (-lim, -dwa options). Note than in case of dead segments (see followings), the automatic procedure might not work and need to be done manually.

Note that for very bad crystals, typically neutron damaged detectors, the default fwhm parameter needs to be increased (the default value is 10) to obtain a reasonable fit result.

Once the calibration seems good the results should be inserted in the configuration file (i.e. add the recal.out coefficients to the 5th column of PreprocessingFilterPSA.conf) so that the processing of the data should be done properly. A python script *colupdate.py* is provided to help with this⁷.

```
$ cp PreprocessingFilterPSA.conf PreprocessingFilterPSA.save # Just in case...
$ python colupdate.py PreprocessingFilterPSA.conf recal.out -c 4 13 -o
PreprocessingFilterPSA.conf
```

- -c column you want to use from file1 and file2 starting from 0 (5th column in PreprocessingFilterPSA.conf (4 starting from 0), 14th column in recal.out (13 starting from 0))
- -o file output

TkT Calibration check

Once all crystals have been calibrated, we need to process a local level replay of the Preprocessing actor to check the calibrations. For this, we need to create a dedicated folder containing a link to the raw data, the Conf folder (containing the new PreprocessingFilterPSA.conf files), the gen_conf.py script adapted to offline femul replay (see section Replay) and a Topology file similar to this:

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 3    CRY
Producer   CrystalProducerATCA
Filter     PreprocessingFilterPSA
Consumer   None
ENDLOOP
```

The “None” consumer do not write any output file (we only need here the spectra files).

For preparing the following steps, we can modify in the gen_conf.py file the WriteDataMask value to 8 to write the “bdat” files used for x-talk corrections. If you don’t plan to perform x-talk corrections, set this value to 0.

The following lines allows to process the replay:

Note: You need to be in the global folder (where Out, Data, Conf are located)

⁷ See zUseful directory in the agapro package

```
$ python gen_conf.py (will generate the replay Conf and Out folders)
$ femul Topology.conf
```

Open the Preprocessing spectrum file:

```
Out/$CryId/Prep__2-40-16384-UI__Ener.spec
```

with TkT to check the energy calibrations (apply the gain factor written in PreprocessingFilter.conf file to obtain the correct energy) .

**** Graphical method with GammaWare ****

The SpectraViewer toolkit of the GammaWare Watchers gives the possibility to open/display the spectra created by the actors in the ROOT environment. The *RecalEnergy* (see previous paragraph) program has been implemented in this toolkit in order to allow for a direct visualization of the energy calibration procedure.

The SpectraViewer toolkit is loaded from the Watchers directory:

```
$ cd /path/to/gammasoftware/LYON/gw/demos/adf
$ root GANILLoadWatchers.C
```

For the energy calibration tool, one needs to load the Ampli.spec producer files of each crystal. In this example, the runs are located in the directory `"/agatadisk/e705/e705"` :

```


analys@gagata-analysis-1-$ ls /agataadocs/e/os/er/
run_0030.dat.23-05-17_17h25n24s
run_0061.dat.23-05-17_17h27m25s
run_0092.dat.23-05-17_17h27m38s
run_1004.dat.23-05-17_01h36n59s
run_1032.dat.30-05-17_20h34m00s
run_0032.dat.23-05-17_17h08h04s
run_0062.dat.23-05-17_21h41n05s
run_0093.dat.23-05-17_05h55m22s
run_1006.dat.23-05-17_02h24n12s
run_1034.dat.31-05-17_02h47m40s
run_0033.dat.21-05-17_19h34m42s
run_0063.dat.23-05-17_04h20m42s
run_0094.dat.23-05-17_07h56m26s
run_1007.dat.23-05-17_05h12m50s
run_1035.dat.31-05-17_05h44m00s
run_0034.dat.21-05-17_20h21n24s
run_0064.dat.23-05-17_04h30n15s
run_0095.dat.23-05-17_08h54m48s
run_1008.dat.23-05-17_08h07h35s
run_1036.dat.31-05-17_08h08h35s
run_0035.dat.21-05-17_22h20m29s
run_0065.dat.23-05-17_09h30m11s
run_0096.dat.23-05-17_11h37m59s
run_1009.dat.23-05-17_11h37m59s
run_1037.dat.30-05-17_11h32m05s
run_0036.dat.22-05-17_08h43h40s
run_0066.dat.23-05-17_08h43h40s
run_0097.dat.23-05-17_10h10m30s
run_1010.dat.23-05-17_12h28n41s
run_1038.dat.31-05-17_12h45n58s
run_0037.dat.22-05-17_12h20m53s
run_0067.dat.23-05-17_11h49n43s
run_0098.dat.23-05-17_11h49n43s
run_1011.dat.23-05-17_11h49n43s
run_1039.dat.31-05-17_12h41n15s
run_0038.dat.22-05-17_12h28h00s
run_0068.dat.23-05-17_12h28h00s
run_0099.dat.23-05-17_13h17m02s
run_1012.dat.23-05-17_17h35n20s
run_1040.dat.31-05-17_17h35n20s
run_0039.dat.23-05-17_13h02m09s
run_0069.dat.23-05-17_13h02m09s
run_0100.dat.23-05-17_13h59n29s
run_1013.dat.23-05-17_13h59n29s
run_1041.dat.31-05-17_13h59n29s
run_0040.dat.23-05-17_12h28n11s
run_0070.dat.23-05-17_12h28n11s
run_0101.dat.23-05-17_16h31n14s
run_1014.dat.23-05-17_16h31n14s
run_1042.dat.31-05-17_21h26m35s
run_0041.dat.23-05-17_06h21n01s
run_0071.dat.23-05-17_06h21n01s
run_0102.dat.23-05-17_16h46m58s
run_1015.dat.23-05-17_02h10m52s
run_1043.dat.01-06-17_03h03m41s
run_0043.dat.23-05-17_06h47n50s
run_0073.dat.23-05-17_13h03n38s
run_0103.dat.23-05-17_19h53m34s
run_1016.dat.23-05-17_04h46n29s
run_1044.dat.01-06-17_03h31m65s
run_0044.dat.23-05-17_06h58m16s
run_0074.dat.23-05-17_13h12n34s
run_0104.dat.23-05-17_20h06m04s
run_1017.dat.23-05-17_04h59m05s
run_1045.dat.01-06-17_03h06m27s
run_0045.dat.23-05-17_07h25m14s
run_0075.dat.23-05-17_13h12n34s
run_0105.dat.23-05-17_20h06m04s
run_1018.dat.23-05-17_05h58m37s
run_1046.dat.01-06-17_12h15m49s
run_0046.dat.23-05-17_07h25m14s
run_0077.dat.23-05-17_13h32n25s
run_0106.dat.23-05-17_22h23h53s
run_1019.dat.23-05-17_05h58m37s
run_1047.dat.01-06-17_12h15m49s
run_0047.dat.23-05-17_07h52m12s
run_0078.dat.23-05-17_13h44n36s
run_0108.dat.23-05-17_01h44n29s
run_1020.dat.23-05-17_05h58m41s
run_1048.dat.01-06-17_12h15s49s
run_0048.dat.23-05-17_08h18h49s
run_0079.dat.23-05-17_13h54m45s
run_0109.dat.23-05-17_01h49m15s
run_1021.dat.23-05-17_07h50m16s
run_1049.dat.01-06-17_13h12n35s
run_0049.dat.23-05-17_08h18h49s
run_0080.dat.23-05-17_13h54m45s
run_0110.dat.23-05-17_02h10m52s
run_1022.dat.23-05-17_07h50m16s
run_1050.dat.01-06-17_13h12n35s
run_0050.dat.23-05-17_09h44n47s
run_0081.dat.23-05-17_14h51n23s
run_0110.dat.23-05-17_05h58m06s
run_1023.dat.23-05-17_14h13n48s
run_1051.dat.01-06-17_14h27m05s
run_0051.dat.23-05-17_10h38n38s
run_0082.dat.23-05-17_10h38n38s
run_0111.dat.23-05-17_09h14n51s
run_1024.dat.23-05-17_17h12n10s
run_1052.dat.01-06-17_15h11n44s
run_0052.dat.23-05-17_11h26n38s
run_0083.dat.23-05-17_15h12n49s
run_0112.dat.23-05-17_12h14n35s
run_1025.dat.23-05-17_18h02m05s
run_1053.dat.01-06-17_15h11n44s
run_0053.dat.23-05-17_11h26n38s
run_0084.dat.23-05-17_15h12n49s
run_0113.dat.23-05-17_12h14n35s
run_1026.dat.23-05-17_18h02m05s
run_1054.dat.01-06-17_15h11n44s
run_0054.dat.23-05-17_13h52h07s
run_0085.dat.23-05-17_15h20n33s
run_0114.dat.23-05-17_16h45n29s
run_1027.dat.23-05-17_22h04n36s
run_1055.dat.01-06-17_14h24n45s
run_0055.dat.23-05-17_14h10n15s
run_0086.dat.23-05-17_16h20n24s
run_0115.dat.23-05-17_17h51n08s
run_1028.dat.23-05-17_04h23n24s
run_0056.dat.23-05-17_14h47h59s
run_0087.dat.23-05-17_18h09n20s
run_0116.dat.23-05-17_23h10m55s
run_1029.dat.30-05-17_02h26m35s
run_0057.dat.23-05-17_17h50m34s
run_0088.dat.23-05-17_18h09n20s
run_0117.dat.23-05-17_23h10m55s
run_1030.dat.30-05-17_02h26m35s
run_0058.dat.23-05-17_16h29n32s
run_0089.dat.23-05-17_20h05n41s
run_1002.dat.23-05-17_00h02n34s
run_1031.dat.30-05-17_18h03h56s
run_0059.dat.23-05-17_17h00m46s
run_0090.dat.23-05-17_20h04n42s
run_1003.dat.23-05-17_01h40n44s
run_1031.dat.30-05-17_18h12n25s
analys@gagata-analysis-1-$

```

This is done using the `Import(TString Include, TString Exclude, TString RunPath)` method:

```
$ Import("Prod*Ampli.spec","","/agatadisk/e705/e705/run_1052.dat.01-06-17_15h11m44s/")
```

The `gen_conf.py` file must then be selected in the pop-up window. It will be used to get the crystal list and the different gain to apply to each spectra.



The ROOT TBrowser will then contain 4 folders :


- Tasks: Used to apply operations on the spectra like energy calibrations
- Tools: Visualization Tools (not used here)
- SpectraViewer: Folder containing all the loaded actor's spectra
- ROOT Files: To get an access to the open ROOT files

To Start the calibration procedure for a given crystal, one first need to plot the crystal spectra (segments and core). For this, double left click on :

Tasks/RunName/LLPWatchers/CrystalName/Prod/H1D_Ampli

Right click on H1D_Ampli:

- ShowSpectra: Plot the crystal map (segments and cores)
- ShowHist: Plot the spectrum of a selected segment (to specify in the pop-up window)
- AxisTools: Some tools to modify the range (and log scale) of the plotted pads.
- GWRecal: Calibration tools (see below)
- WriteCalibFiles: to store the calibration output in the PreprocessingFilterPSA.conf files
- PropagateTasks: not used here
- ShowClassInfo: not used here
- SetLoupe: if boolean set to true in the pop-up window, allow to zoom on a selected pad



For the energy calibration, the second library (lib 1 starting from 0) is used (see explanations in the previous part). Use the ShowSpectra method to plot all the spectra in the crystal map. The axis tool, FullRangePerPad (AxisTools/Range/FullRangePerPad) can be used to define an individual automatic range to each pad.

To specify the calibration source, use :

GWRecal/Configure/TheoreticalPeaks/SetSource


Defined sources: 22Na 40K 56Co 57Co 60Co 88Y 133Ba 134Cs 137Cs 152Eu 208Pb 226Ra 241Am (more than one source can be set, separated by a space).

The keyboard shortcut “s+s” can also be used (for Set Source)

To fit the whole crystal map with the default parameters, use:

```
GWRecal/Calibrate/FitAll ("Ctrl+f" shortcut)
```

The SetLoupe method can here be useful to carefully check the fit result on the different segments, when active, select a pad by a wheel click and then press the space key.



A window range can be set for the peak search using :

```
GWRecal/Configure/Channel/SetGlobalChannelLimits
```

```
"Ctrl pressed and c+l"
```

=> For all pads

```
GWRecal/Configure/Channel/SetHistChannelLimit
```

```
"c+l"
```

=> To be applied only on a specified segment (in this case, the shortcut applies for the selected pad (mouse wheel)).

Peak limits (minimum FWHM and maximum) can also be defined

```
GWRecal/Configure/PeakSearch/SetGlobalPeaksLimits
```

```
"Ctrl pressed and p+l"
```

=> For all pads

GWRecal/Configure/PeakSearch/SetHistPeaksLimits

"p+1"

=> To be applied only on a specified segment (in this case, the shortcut applies for the selected pad (mouse wheel)).

When new parameters have been set for the peak search, the whole crystal map can be re-calibrated ("Ctrl+F"), or only the selected pad ("F"). Once all segments have been fine tuned, make a final fit of the whole map to print the final fit parameters. Fit results are also written on disk in the FitResults/GWRecal directory.


```
# Spectra taken from file : /media/dudouet/ExternalDisk/Experiments/e705/Eu_Run_0036/Data/00A/Prod_4-38-32768-UI_Ampli.spec
# lib #spec #pks #ok rEnergy FW05 FW01 Area Position Width Ampli WTML WTMR slope*gain rChi2%
#
1 0 28 10 1408.20 3.723 9.316 2042 9329.24 20.4 66 4.221 1.823 0.301890 3.86
1 1 19 10 1408.93 4.572 13.721 1009 8975.26 12.5 25 11.978 1.990 0.313958 18.84
1 2 19 10 1408.39 4.933 13.459 1422 9634.25 22.2 32 6.471 1.823 0.292371 7.09
1 3 19 10 1408.38 3.980 10.808 1358 8803.02 16.6 42 6.298 1.823 0.319975 3.75
1 4 10 10 1408.27 4.085 10.717 778 8739.85 18.8 24 5.239 1.823 0.322263 8.70
1 5 8 7 1408.41 3.185 7.905 384 8608.14 16.4 16 4.080 1.823 0.327229 9.78
1 6 24 10 1408.32 3.669 8.979 1919 9340.70 21.0 64 3.860 1.823 0.301546 2.37
1 7 19 10 1408.38 4.049 7.381 725 8810.94 25.3 27 1.823 1.823 0.319688 10.41
1 8 20 10 1409.18 4.530 13.186 1375 9337.56 14.6 34 10.145 1.823 0.301831 31.17
1 9 12 10 1407.71 4.678 8.526 1049 9416.49 31.3 31 1.823 1.823 0.298988 4.28
1 10 12 8 1408.08 4.157 7.577 593 9127.03 26.9 21 1.823 1.823 0.308553 13.62
1 11 8 7 1408.35 3.232 7.620 372 9048.41 18.9 15 3.366 1.823 0.311292 8.45
1 12 29 10 1408.46 3.698 9.453 1753 8517.95 17.7 62 4.630 1.823 0.330703 6.42
1 13 20 10 1408.98 4.461 12.172 874 8878.49 18.5 24 6.471 1.823 0.317392 30.38
1 14 18 10 1409.01 5.976 17.490 1691 9438.52 18.8 31 10.628 1.823 0.298566 25.96
1 15 18 10 1407.95 5.492 10.014 1055 9145.88 35.7 28 1.823 1.824 0.307887 15.39
1 16 14 10 1407.94 3.621 6.599 512 9430.52 24.3 20 1.823 1.823 0.298593 9.98
1 17 11 8 1408.39 2.533 6.187 329 9084.28 14.1 16 3.830 1.823 0.310072 8.02
1 18 22 10 1408.37 4.024 10.337 1231 9039.84 20.2 38 4.738 1.823 0.311592 10.29
1 19 18 10 1407.93 5.344 9.741 512 8941.51 33.9 14 1.823 1.823 0.314921 20.93
1 20 18 10 1408.52 5.377 14.457 939 9217.08 24.3 21 5.951 1.823 0.305632 20.48
1 21 17 10 1407.31 4.748 8.653 626 8477.21 28.6 21 1.823 1.823 0.332022 11.43
1 22 11 10 1408.22 3.162 8.002 416 9179.04 16.7 16 4.430 1.823 0.306833 4.20
1 23 8 7 1407.90 2.483 5.322 193 9441.83 16.4 10 2.536 1.823 0.298225 7.77
1 24 29 10 1408.47 3.684 9.510 1223 9197.60 18.6 40 4.842 1.823 0.306270 10.56
1 25 19 10 1407.50 5.207 9.564 596 8675.48 32.1 17 1.823 1.851 0.324478 15.92
1 26 20 10 1408.32 4.789 12.565 1120 8373.87 21.2 31 5.239 1.823 0.336362 28.35
1 27 18 10 1407.93 4.754 12.921 1291 9263.10 20.9 32 6.322 1.823 0.303988 9.36
1 28 15 10 1408.05 3.318 7.994 683 9161.08 19.1 26 3.634 1.823 0.307398 3.33
1 29 13 8 1407.92 3.463 7.435 391 9628.30 23.3 14 2.546 1.823 0.292454 1.88
1 30 31 10 1408.19 3.724 9.666 1768 9867.09 19.9 53 4.971 1.823 0.285431 2.71
1 31 17 10 1408.19 5.282 14.543 1171 9186.03 21.9 26 6.830 1.823 0.306594 26.38
1 32 17 10 1407.79 4.516 8.230 1099 9106.31 29.2 35 1.823 1.823 0.309191 4.54
1 33 15 10 1407.74 4.470 8.147 1039 9751.79 31.0 32 1.823 1.823 0.288714 6.03
1 34 15 9 1408.20 3.668 9.377 813 9244.93 19.1 27 4.633 1.823 0.304643 9.90
1 35 12 7 1407.41 3.630 6.617 393 9217.33 23.8 16 1.823 1.823 0.305383 10.92
1 36 52 10 1407.99 3.296 6.772 40858 5931.57 13.8 2590 2.298 1.823 0.474745 1.11
1 37 23 10 1408.06 3.481 7.502 40894 1652.52 4.0 8602 2.572 1.823 1.704130 274.27
Calibration File write in: FitResults/GWRecal/CurrentRun/00A_Amplitude_Lib1.fit
```

For any reminder on the possible commands and keyboard shortcuts:

GWRecal/PrintHelp ("p+h")

Once a crystal map has been correctly fitted, a "PlotResults" tool allows to print the different parameters versus the segment ID, or versus other parameters :

GWRecal/PlotResults/PlotCalibResults ("p+c")



A detailed explanation of this PlotResult tool can be shown using :

```
GWRecal/PlotResults/PrintPlotHelp
```

Finally, a global fit of the whole crystals can be done in a single command, using the PropagateTask method on the “LLPWatchers” folder. A detailed help on this tool can be displayed by a right click on LLPWatchers folder and :

```
PropagateTask/PrintHelp
```

The PropagateTask tool can be used using :

```
PropagateTask/PropagateTask
```

For example, the command :

```
"ShowSpectra Lib=1 ; Recal Source=152Eu PL=20,10 ; Plot  
Exp=FW05:spec Opt=same ; all"
```

will plot the first library of all segments of all crystals. Spectra will then be fitted with the Recal code using an ^{152}Eu source and limiting the peak search to peaks FWHM < 20 and amplitude > 20. FWHM of the reference peak for each crystal will finally be plotted on a same Canvas.

To store the fit results in the PreprocessingFilterPSA.conf file, you can use the colupdate.py script (see the previous section) using the files written in the FitResults folder, or use the WriteCalibFiles option (right click on a crystal), and giving the path of the Conf folder (it will be applied for the whole crystals).

Pole-zero / shaping-time adjustments

Explanations and goals

The **trapezoidal filter** is applied in the pre processing actor. We have here the possibility to change the shaping time and tune the decay time of the exponent for the experiment (depending on the rate/crystal).

To distinguish the meaning of pole-zero and shaping-time, in the calibration scripts:

- **pole-zero** is related to the **decay-time**
- **shaping-time** to the **rise-time**

In order to set the pole-zero of the cores (for the segments one should leave them like they are) it is necessary to fit **the exponential of the long traces** (this is done on site by the local team when a new crystal is added in the topology).

NB: this part is done on site at the beginning of an experiment and should not be modified afterwards

Access to the following account-machine from the AGATA server:

```
$ ssh psa-tests@scgw3
$ cat config_crystal/crystal_${CryId}.rc
```

```
...
all risetime 500
all baseline 4
all polezero 4600
CC channel 1
CC threshold 100
CC polezero 4606
...
```

These coefficients for pole zero and rise should be integrated in the PreprocessingFilterPSA.conf.

Crosstalk corrections

Explanations and goals

Crosstalk is present in any segmented detectors. Through couplings, the collection of the signal in one segment modifies the signal collected in the neighboring ones. One consequence is a shift in energy (applying a calibration based on events with just one segment fired) which increases as a function of the number of segment fired in one crystal. There are different crosstalk in segmented Germanium crystals, in particular differential and proportional. To have more information on such effect in different segmented Germanium detectors, you can read the following papers (non-exhaustive list):

Vetter NIMA 452 (2000), Swensson NIMA 540 (2005), B.Rossé NIMA 565 (2006), B. Bruyneel NIMA 608 (2009)

Concerning the AGATA capsule, it is possible use the crosstalk correction procedure to recover up to one broken or missing segment per crystal (by using the fact that sum of the energies of the segments should be equal to the core energy).

Tools available

**** TkT and command line programs ****

In order to perform this correction, we will use the files with the extension *.bdat (**event_energy.bdat**) produced by the CrystalProducer if the 4th bin in the WriteDataMask instruction of *CrystalProducer.conf* is set ("WriteDataMask 8"). The program used to perform this is called **SortCapsule**. Xtalk corrections needs to be done using a ^{60}Co source.

1) First of all, it is necessary to write the preprocessing calibration file in a new format, use the following command to prepare the file in the proper format for the next steps:

```
$ cd Conf/${CryId}          # Go in the working directory if not already done
$ rm -f ecalF1.cal          # if file was already existing

$ for i in {1..38} ; do cat PreprocessingFilterPSA.conf | head -n $i | awk -v var=$(( $i - 1 )) '{print "0\t" var "\t2\t0\t" $5}' | tail -n 1 >> ecalF1.cal ; done ;
```

The ecalF1.cal file should be similar to:

0	0	2	0	0.301769
0	1	2	0	0.313686
...				
0	36	2	0	0.484332
0	37	2	0	1.740927

2) Verify the calibration coefficients

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -egain 5
```

this will generate the following files for each crystal:

```
- proj__3-38-32768-UI__raw.spec raw projections      : [0]=Tot [1]=SG@F1 [2]=CC@F1
- spec__3-38-16384-UI__cal.spec calibrated spectra  : [0]=Tot [1]=SG@F1 [2]=CC@F1
- ssum__2-12-16384-UI__cal.spec sumEnergy[segFold]  : [0]=SG [1]=CC
```

with [N] the library number.

The two last files can be read using TkT to check the segments energies. Energies needs to be good for fold 1, but note that the position of the peaks in the first set of spectra in ssum__2-12-16384-UI__cal.spec depends on the number of fired segments (Fold). No such dependency for the second set (core).

3) Call xTalkSort to sort and analyze the AGATA events dumped into event_energy.bdat.0000 by producing the cross-talk spectra for full energy release of 1332.5 keV in one segment (F1)

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -egain 5
-specXT -trigewin 1325 1340
```

With :

- ecalF1: file with energy calibration coefficients for singles extracted above
- egain: gain factor for energy spectra
- specXT: generate cross talk spectra from Mseg=1
- trigewin: energy window on trigger channel

4) Once the different spectra have been created, the RecalEnergy program can be used to build the direct Xtalk matrix:

```
$ RecalEnergy -spe xspe__36-37-16384-UI__cal.spec -num 1332 -ener 1332.5 -gain 5 -offs
1000 -Xtalk 37 |tee 1325-1340.txt
```

With :

- Xtalk 37: Use the 37th spectrum for the core in the Xtalk calculation
- offs 1000: channel offset to subtract to the position of the peaks

5) Transform this file to the proper format for direct cross talk coefficients:

```
$ grep -v "^#" 1325-1340.txt |grep -v "^ *36 " |cut -b15-102 --complement |tee
xdir_1325-1340.cal
```

6) Invert the cross talk coefficients and rename it:

```
$ xTalkInvert -f xdir_1325-1340.cal ; mv xdir_1325-1340.cal.inv xinv_1325-1340.cal
```

The inverse matrix is given the name of xdir_1325-1340.cal.inv. The traditional name of this file was called xinv_1325-1350.cal. Up to you to rename the file or keep the given name.

7) Call again xTalkSort using the Xtalk inverted matrix to check the corrected spectra:

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -egain 5
-recalXT xinv_1325-1340.cal
```

This will produce the spectra files:

```
- proj__3-38-32768-UI__raw.spec raw projections      : [0]=Tot [1]=SG@F1 [2]=CC@F1
- spec__3-38-16384-UI__cal.spec calibrated spectra  : [0]=Tot [1]=SG@F1 [2]=CC@F1
- spec__3-38-16384-UI__adj.spec adjusted spectra    : [0]=Tot [1]=SG@F1 [2]=CC@F1
- ssum__2-12-16384-UI__cal.spec sumEnergy[segFold]  : [0]=SG [1]=CC
- ssum__2-12-16384-UI__adj.spec sumAdjusted[segFold]: [0]=SG [1]=CC
```

with [N] the library number.

Check with TkT the spectra in spec__3-38-16384-UI__adj.spec and ssum__2-12-16384-UI__adj.spec. In particular the shift in position of the sum energy of the segments as a function of the segment fold should have disappeared.

Specific case of “lost” or “broken” segments:

The part extends the treatment to the case of one **dead** segment which can be “lost” or “broken”.

A **broken** segment is the result of a problem at the FET level with the consequence that the charge of the segment is not collected but flows to the neighbors. The salient effect of a **broken** segment is the presence of ghost peaks in the neighbors and of a strong step-like tail in the spectrum of the core. The ghost peaks and the left step can be seen as enhanced cross-talk.

The case of a **lost** segment is when the detector works normally but the signal is not present in the data due, e.g., to a broken wire or a faulty digitizer channel. In this case there are no ghost peaks.

It is worth remarking that segments with **unstable** gain could be transformed into (and treated as) **lost** segments by setting their energy calibration to zero. Of course this is possible only if all other segments in the detector work correctly.

To treat **dead** segments, we have to find a way to quantify the amount of missing energy. For **broken** segments, we have also to generate a specific set of cross-talk correction coefficients capable of removing the ghost peaks from the affected neighbors. Two different ways of treating this are possible

1) In the calibration file adapted for Xtalk calculations (ecalF1.cal) **and** in the PreprocessignFilterPSA.conf, the calibration parameter of the **dead** segment **must to be set to 0**.


(For the command lines in the following, it is supposed that you are in the Conf folder of the concerned detector named in the following: "\$CryId")

2) To determine the energy released in the missing segment:

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -matCCSG
```

Open with Mat.exe the file CC-SG__50-1500-1500-US__ma.matr. Matrix #36 shows the correlation between the energy seen by the core and the sum-energy of all segments.

For a **broken** segment it should look like the following:



The main diagonal is not at 45° (as usual when the sum-energy is over all segments) and there is a lower limit line for the lost in the sum of segments (because part of the energy released in the **broken** segment is “collected” in its close neighbors, generating the “ghost peaks”) and also a small loss of core energy (the slightly tilted lines coming down from the peaks which would be perfectly vertical if the segment would be disconnected).

The coefficients to manage the broken segment (**here and also in PreprocessingFilterPSA**) are:

- **deadXsg** = $1267.5/1332.5 = 0.95122$
→ **slope of main diagonal**
- **deadXcc** = $(1332.5-1312.9)/(1267.5-366.2) = 0.021746$
→ **1/slope of core loss**

For a **lost** segment:

the main diagonal has practically the same slope but the lower line is not present and the lines out of the peaks are just vertical. The coefficient to manage the lost segment are:

- **deadXsg** = 0.95122
→ **slope of main diagonal**
- **deadXcc** = 0
→ **no core loss**

3) We can verify these values from the core-segment correlation matrices produced after giving the index of the dead segment and the two coefficients:

For a **lost** segment (case of dead segment id = 4 [A5]):

We will use in the following the environment variable \$DeadSeg=4

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -deadSeg  
$DeadSeg 0.95122 0. -matCCSG
```

The off diagonal stuff in sCC-eSumSG has disappeared by construction because the missing energy has been assigned to the lost segment. Indeed A5 is now present with the correct energy and we are essentially back to the normal situation and we can produce the cross-talk matrix as:

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -deadSeg  
$DeadSeg 0.95122 0. -matx1
```

and the cross talk correction coefficients as in the normal case:

```
$ xTalkMake -f xSG__36-36-100-1536-US__ij.matr
```

Test the result (see step 7 of standard Xtalk corrections):

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -deadSeg  
$DeadSeg 0.95122 0. -recalXT xinv_1325-1340.cal
```

For a broken segment (case of dead segment id = 4 [A5]):

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -deadSeg  
$DeadSeg 0.95122 0.021746 -matCCSG
```

Also for this case the off diagonal stuff in sCC-eSumSG has disappeared and the matrix of the broken segment is no more empty. However, the assigned energy is much too small due to the fact that part of it is still in the neighbors. The slope of the diagonal line in the two matrices is:

- **slopeDeadSG** = $897.5/1332.5 = 0.6734$

→ (matrix #4) **needed for xTalkMake**

To correct the remaining effects we need to determine the proper xTalk coefficients. The procedure described in the first part of this document is not suited as it looks only for slightly negative-energy peaks, missing completely the positive-energy ghosts. The program manages this case from the fact that the third parameter of the command line switch “-deadSeg” is positive:

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -egain 5  
-deadSeg $DeadSeg 0.95122 0.021746 -matx1
```

In the generated file : xSG__36-36-1000-100-US__ij.matr, like for the standard xTalk matrices the x axis (second-last index) is the spectrum of the “affected segment” while the y axis (last index) is the spectrum of the “affecting (net-charge) segment”. In this case the 100 channels-long Y axis is a [1300-1350] keV window on the core energy (0.5 kev/chan), while the 1000 channels-long x-axis is the spectrum of the affected segment for the energy range [-15 185]-keV (at a gain of 5

channels/keV). The matrix is incremented only if one segment sees more than 60% (>800 keV) of the core energy and that segment is considered as the affecting one. Have a look to the matrices [4] [seg] to see the location of the ghost peaks in close neighbours of A5. Then, calculate the cross-talk correction coefficients by modifying the diagonal element of the broken segment with:

```
$ xTalkMake -f xSG__36-36-1000-100-US__ij.matr -dxy 1000 100 -gate 55 75 -xcal -15 5
-egam 1332.5 -dval $DeadSeg 0.6734
```

To be consistent with the Xtalk file names, rename the matrices:

```
$ mv xdir_0055-0075.cal xdir_1325-1340.cal
$ mv xinv_0055-0075.cal xinv_1325-1340.cal
```

Finally, test the result with (see step 7 of standard Xtalk corrections):

```
$ xTalkSort -ifile ../../Out/$CryId/event_energy.bdat.0000 -ecalF1 ecalF1.cal -deadSeg
$DeadSeg 0.95122 0.021746 -recalXT xinv_1325-1340.cal
```

4) In the *gen_conf.py*, define in the PreprocessingFilter and PSA actors the **dead** segments as in the example bellow with 00A as a broken segment, and 11A as a lost one:

```
#####
PreprocessingFilter=(
  "ActualClass      PreprocessingFilterPSA", # name of the used daughter class
  "SaveDataDir      $SAVEDIR/$CRYSTAL",      # normally Out/01A...
  "EnergyGain       4",                      # channels/keV of the calibrated energy spectra
  "XtalkFile        xinv_1325-1340.cal",      # cross talk correction coeffs for the energies
  "WriteTraces      100",                   # number of traces written
  ##### command lines to be produced only for the specified crystals
  {
    '00A' : ("DeadSegment    29 0.950 0.049"), # case of a broken segment
    '11A' : ("DeadSegment    18 0.947 0."),    # case of a lost segment
  }
)
#####
PSAFilter=(
  "ActualClass      PSAFilterGridSearch",    # name of the used daughter class
  "BasisFile        $SIGNAL_BASIS",          # this is generated from the GeDataBase structure
  "SaveDataDir      $SAVEDIR/$CRYSTAL",      # normally Out/Data(online)
  "EnergyGain       4",                      # channels/keV of the calibrated energy spectra
  "XtalkFile        xdir_1325-1340.cal",      # cross talk correction coeffs for traces
  "Threads          5 300",                  # number of threads, events/thread
  "GridSearchType   Adaptive",               # SegCenter, Adaptive, CoarseOnly or Full;
  ##### command lines to be produced only for the specified crystals
  {
    '00A' : ("DeadSegment    29"),
    '11A' : ("DeadSegment    18"),
  }
)
)
```

Check also that the files: xinv_1325-1340.cal xdir_1325-1340.cal are written in the “ExtraFiles” list of the *gen_conf.py*:

```
ExtraFiles={
  'CRYSTAL' : "CrystalProducerATCA.conf PreprocessingFilterPSA.conf xinv_1325-1340.cal xdir_1325-1340.cal",
}
```

Time alignment

Explanations and goals

The Preprocessing has another very important “ingredient”, i.e. the possibility to aligned in time (based on Straight line fit or Digital CFD of the signal rise-time) the traces of the different segments to the core one.

This operation is very important to obtain good performances of the PSA algorithm. In other words, the PSA tends to ‘process’ all the waveforms as beginning at the same time and therefore, if there is any mis-match in time, the reconstructed interaction points are prone to be unreliable. A bad time alignment can result in clusterisation patterns at the output of the PSA.

Tools available

**** command line programs ****

The spectra to be used are in the *Prep__6-40-1000-UI__TT.spec*⁸ (before was *Prep__2-40-1000-UI__TT1.spec*). Library 1 (starting from 0) corresponds to timing before applying any shift. It must be used for time calibration. Library 3, after time shift, can be used to check the quality of the time alignment.

1) Time alignment of segments:

The command lines to be used per crystal is the following:

```
$ cd Conf/$CryId # If not already done, go in the working directory
$ RecalEnergy -spe ../../Out/$CryId/Prep__6-40-1000-UI__TT.spec -sub 40 -num 36 -T 500 |
tee shift_TT.out
```

Options:

- -sub 40: Second library
- -num 36: 36 segments
- -T 500: we want all the peaks at this position

Add this coefficients to the 7th column of *PreprocessingFilterPSA.conf* using the *colupdate.py* script:

```
$ cp PreprocessingFilterPSA.conf PreprocessingFilterPSA.save2 # Just in case...
$ python colupdate.py PreprocessingFilterPSA.conf shift_TT.out -c 6 13 -o
PreprocessingFilterPSA.conf
```

2) Time alignment of the core:

The time alignment of the core is done using the spectrum file produced by the PSA actor name “*Psa__40-1000-UI__Tzero.spec*” (Here we don’t care about what the PSA actor is doing, we only use this spectrum file). **A replay from traces to PSA is thus required, taking into account the new time calibration of the segments.**

Note: Be careful to have in the *gen_conf.py* file, the option “NoMultiHist” **commented** in the PSA actor (otherwise the spectra files will no be produce).


⁸ This set of spectra is produced by the Preprocessing filter.

This kind of topology can be used:

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 4    CRY
Producer   CrystalProducerATCA
Filter      PreprocessingFilterPSA
Filter      PSAFilter
Consumer    None
ENDLOOP
```

The statistic needed here is not very important, few minutes should be enough. To align the timing of the core, look at the histogram n° 36 of the file Psa__40-1000-UI__Tzero.spec as shown below:



A well align core must contain only one peak, centered at channel 500. In the example above, a broad peak is observed with a peak position at 478. A shift of 22 channels is thus needed. Apply this shift to the timing value of the two cores set in the PreprocessingFilterPSA.conf file:


Before core alignment:

```
...
core          0  4390  500  0.482203  0  5.
core          1  4390  500  1.735468  0  5.
tntf  2097152
```

After core alignment:

```
...
core          0  4390  500  0.482203  0  27.
core          1  4390  500  1.735468  0  27.
tntf  2097152
```

Once the core has been aligned, start a new short replay to check the result, it should be similar to:



Finally, reprocess a replay to check the quality of the time alignment, looking at the spectra file “Prep__6-40-1000-UI__TT.spec” (library 3)

The Pulse Shape Analysis (PSA) actor

The PSA actor is used to extract from the signals shapes the position of each interaction point with a 5mm precision. No specific calibration are needed for this step. But the quality of the PSA is highly dependent on the good calibrations at the Preprocessing level.

The PostPSA actor

The PostPSA filter actor allows to make all the final operations on the local level data (neutron damage correction, final energy calibration with an offset and core time alignment).

neutron damage corrections

The damages caused by interactions between neutrons and the Germanium detectors deteriorate the gamma spectra quality. The typical effect of this neutron damages is a left tail on the peaks. This detector deterioration is increasing along time and its correction is mandatory to obtain a satisfying energy resolution.

This corrections is done in three steps:

1. First, a recalibration of the segments and core can be performed before the neutron damage correction.
2. The program named “SortPsaHits” is then used to estimate the neutron damages and correct the energies.
3. After the neutron damage correction, a final recalibration is processed to correct from possible shifts induced by the neutron damage correction.

0) First, it is necessary to do a replay of a ⁶⁰Co source run from the traces including the

“WritePsaHits” option in the PSA actor of the *gen_conf.py* file:

Note: It is necessary to have configured the correct path to the ADL signal bases in the *gen_conf.py* (see the Replay section)

```
#####
PSAFilter=(
"ActualClass      PSAFilterGridSearch",      # name of the used daughter class
"BasisFile        $SIGNAL_BASIS",            # this is generated from the GeDataBase structure
"SaveDataDir      $SAVEDIR/$CRYSTAL",        # normally Out/Data(online)
"EnergyGain       4",                        # channels/keV of the calibrated energy spectra
"XtalkFile        xdir_1325-1340.cal",        # cross talk correction coeffs for traces
"Threads          5 300",                    # number of threads, events/thread
"GridSearchType   Adaptive" ,               # SegCenter, Adaptive, CoarseOnly or Full;
"WritePsaHits" ,                               # writes the hits in binary
)
```

This will produce the spectra files: Psa__0-16-F__Hits.fdat, needed to determine the coefficients of the neutron damage corrections.

This kind of topology can be used:

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 4   CRY
Producer  CrystalProducerATCA
Filter    PreprocessingFilterPSA
Filter    PSAFilter
Consumer  BasicAFC
ENDLOOP
```

The use of the BasicAFC consumer here will write the psa*.adf files that will be used for faster replays in the following.

In your calibration folder, create a new symbolic link “Out” linked to the output of this new replay

1) Prepare a dummy trapping file:

```
$ cd Conf/$CryId          # Go in the working directory if not already done
$ rm -f Trapping.cal      # If already existing
$ for i in {0..35} ; do echo -e "$i\t1.\t1.\t999999.9\t999999.9\t1.\t1." >> Trapping.cal
; done
```

This file should be like following:

0	1.	1.	999999.9	999999.9	1.	1.
1	1.	1.	999999.9	999999.9	1.	1.
...						
...						
34	1.	1.	999999.9	999999.9	1.	1.
35	1.	1.	999999.9	999999.9	1.	1.

It can be understood as follow:

column 0: segment id

column 1: extra gain on the segment before correction

column 2: extra gain on the core before correction

column 3: electron-trapping correction (lambdaE)

column 4: hole-trapping correction (lambdaH)

column 5: extra gain on the segment after correction

column 6: extra gain on the core after correction

2) The program *SortPSAHits* is necessary to sort PSA hits in a specific format to determine neutron damage correction parameters.

```
$ SortPsaHits -f ../../Out/$CryId/Psa__0-16-F__Hits.fdat -best 1300 1350 -bpar 1 10000 0  
|tee sort_hit.log
```

3) Update the dummy trapping file:

```
$ tail -n 39 sort_hit.log |head -n 36 > sort_hit_nohead.log;  
$ awk 'FNR==NR{a[NR]=$4;next} {$4=a[FNR]}1' sort_hit_nohead.log Trapping.cal |awk  
'{printf "%2s %10s %10s %10s %10s %10s %10s \n", $1,$2,$3,$4,$5,$6,$7}' >  
Trapping_tmp.cal;  
$ awk 'FNR==NR{a[NR]=$5;next} {$5=a[FNR]}1' sort_hit_nohead.log Trapping_tmp.cal |awk  
'{printf "%2s %10s %10s %10s %10s %10s %10s \n", $1,$2,$3,$4,$5,$6,$7}' |tee  
Trapping.cal;
```

The columns corresponding to neutron damage correction should now be updated in the Trapping.cal files like following:

0	1.	1.	51.6	13.4	1.	1.
1	1.	1.	1482.1	11.1	1.	1.
...						
...						
34	1.	1.	326.8	11.1	1.	1.
35	1.	1.	385.5	11.1	1.	1.

4) Generation of the file Pso__2-4-40-2048-UI__Ener.spec for the recalibration of the segment before neutron correction:

```
$ SortPsaHits -f ../../Out/$CryId/Psa__0-16-F__Hits.fdat -gain 5 -offs 5000 -fcal  
Trapping.cal
```

This Pso__2-4-40-2048-UI__Ener.spec file contains two spectra libraries of 4 sub libraries of 40 spectra:

- lib 0 – Segments, lib 1 Core:
 - sublib 0: – original spectra
 - sublib 1: – original spectra + recal before neutron damage correction
 - sublib 2: – spectra after neutron damage correction
 - sublib 3: – neutron correction + final recalibration

5) If necessary, recalibration of the segments and cores before neutron correction (the search peaks parameters “-dwa” probably needs to be adapted) :

- segments recalibration:

```
$ RecalEnergy -spe Pso__2-4-40-2048-UI__Ener.spec -num 36 -sub 0 -gain 5 -offs -5000  
-noTR -dwa 30 5 | tail -n 36 |tee log_sg_pre.cal
```

- cores recalibration:

```
$ RecalEnergy -spe Pso__2-4-40-2048-UI__Ener.spec -num 36 -sub 160 -gain 5 -offs -5000
-noTR -dwa 30 5 | tail -n 36 |tee log_cc_pre.cal
```

6) Insert the pre-trapping recalibration parameters in the trapping file:

```
$ awk 'FNR==NR{a[NR]=$14;next} {$2=a[FNR]}1' log_sg_pre.cal Trapping.cal |awk '{printf
"%2s %10s %10s %10s %10s %10s %10s \n", $1,$2,$3,$4,$5,$6,$7}' |tee Trapping_tmp.cal
```

```
$ awk 'FNR==NR{a[NR]=$14;next} {$3=a[FNR]}1' log_cc_pre.cal Trapping_tmp.cal |awk
'{printf "%2s %10s %10s %10s %10s %10s %10s \n", $1,$2,$3,$4,$5,$6,$7}' |tee
Trapping.cal
```

The columns corresponding to pre-calibration should now be updated in the Trapping.cal files like following:

0	0.999311	1.000675	51.6	13.4	1.	1.
1	1.000594	1.000559	1482.1	11.1	1.	1.
...						
...						
34	0.999558	1.000857	326.8	11.1	1.	1.
35	0.999976	1.000779	385.5	11.1	1.	1.

7) Apply the recalibration:

```
$ SortPsaHits -f ../../Out/$CryId/Psa__0-16-F__Hits.fdat -gain 5 -offs 5000 -fcal
Trapping.cal
```

This will re-generate the file Pso__2-4-40-2048-UI__Ener.spec

8) Post trapping recalibration of the segments and cores after neutron correction (the search peaks parameters “-dwa” probably needs to be adapted) :

- segments recalibration:

```
$ RecalEnergy -spe Pso__2-4-40-2048-UI__Ener.spec -num 36 -sub 80 -gain 5 -offs -5000
-dwa 30 5 | tail -n 36 |tee log_sg_post.cal
```

- cores recalibration:

```
$ RecalEnergy -spe Pso__2-4-40-2048-UI__Ener.spec -num 36 -sub 240 -gain 5 -offs -5000
-dwa 30 5 | tail -n 36 |tee log_cc_post.cal
```

9) Insert the post-trapping recalibration parameters in the trapping file:

```
$ awk 'FNR==NR{a[NR]=$14;next} {$6=a[FNR]}1' log_sg_post.cal Trapping.cal |awk '{printf
"%2s %10s %10s %10s %10s %10s %10s \n", $1,$2,$3,$4,$5,$6,$7}' |tee Trapping_tmp.cal
```

```
$ awk 'FNR==NR{a[NR]=$14;next} {$7=a[FNR]}1' log_cc_post.cal Trapping_tmp.cal |awk
'{printf "%2s %10s %10s %10s %10s %10s %10s \n", $1,$2,$3,$4,$5,$6,$7}' |tee
Trapping.cal
```

The columns corresponding to post-calibration should now be updated in the Trapping.cal files like following:

0	0.999311	1.000675	51.6	13.4	0.998467	0.998754
1	1.000594	1.000559	1482.1	11.1	0.999063	0.999600
...						
...						
34	0.999558	1.000857	326.8	11.1	0.999425	0.999294
35	0.999976	1.000779	385.5	11.1	0.999054	0.999287

10) Re-generation of the file Pso__2-4-40-2048-UI__Ener.spec for final checks (be careful, a gain of 5 and an offset of 5000 is applied, the 1333 peak will thus be located at channel 1665):

```
$ SortPsaHits -f ../../Out/$CryId/Psa__0-16-F__Hits.fdat -gain 5 -offs 5000 -fcal
Trapping.cal
```

Finally, the *gen_conf.py* needs to be adapted in the PostPSA part, as follow (example for crystal 00A):

Note: If this line is already present, with other keywords, **remove them**, they will be added later and **must not be present for the moment**.

```
#####
PostPSAFilter=(
"ActualClass      PostPSAFilter",          # name of the used daughter class
"SaveDataDir      $SAVEDIR/$CRYSTAL",      # normally Out/Data(online)
"EnergyGain       4",                      # channels/keV of the calibrated energy spectra
### command lines to be produced only for the specified crystals
{
'00A' : ("TrappingFile Trapping.cal"),
}
)
```

Check also that the file: Trapping.cal is written in the “ExtraFiles” list of the *gen_conf.py*:

```
ExtraFiles={
'CRYSTAL' : "CrystalProducerATCA.conf PreprocessingFilterPSA.conf xinv_1325-1340.cal xdir_1325-
1340.cal Trapping.cal",
}
```

Final energy re-calibrations (with offsets)

In all the previous steps, energy calibrations were always linear. But the PostPSA actor allows to apply a final segment and energy recalibration, including an offset and if necessary higher order of calibrations. This step needs thus to be done on data taken with an ¹⁵²Eu source to have a good estimation of the non-linearity of the detectors.

0) The first step is to perform a replay of an ¹⁵²Eu run, from the traces, including the PostPsaFilter actor with neutron damage correction applied, and **without recalibration and ForceSegmentToCore** options.

For an easier and faster work in the following, this replay can be done in two steps:

- a) Replay from traces to PSA (without the post PSA) and store the psa*.adf files using the BasicAFC consumer:

You need to update the symbolic link “Data” in order to link to the ^{152}Eu data folder.

```
$ mv Out Out_save          # If you want to save the last replay from  $^{60}\text{Co}$ 
$ rm Data
$ ln -s /path/to/ $^{152}\text{Eu}$ /run/Data Data
$ python gen_conf.py
```

This kind of topology can be used:

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 5    CRY
Producer   CrystalProducerATCA
Filter     PreprocessingFilterPSA
Filter     PSAFilter
Filter     PostPSAFilter
Consumer   BasicAFC
ENDLOOP
```

A large statistics is better, this replay can be long, but it will allow to no more perform the PSA later. The next replays will then be very fast.

- b) Replay of the PostPSA only, using as input the psa*.adf files that have been produced at step a). For this, you need to move the Out folder in which the psa*.adf have been written in a new folder (“Data_Local” for example), and update the symbolic link “Data” to this folder, in such a way that the input files of the replay will be the psa*.adf that have been just produced:

```
$ mv Out Data_Local
$ rm Data
$ ln -s Data_Local Data
$ python gen_conf.py
```

This kind of topology can be used:

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 3    CRY
Producer   BasicAFP
Filter     PostPSAFilter
Consumer   None
ENDLOOP
```

This will produce the “Post__5-40-16384-UI__Ener.spec” spectra filename

- 1) Segments recalibration:

- case 1 : for detectors with Trapping.cal files **present**:

```
$ RecalEnergy -spe ../../Out/$CryId/Post__5-40-16384-UI__Ener.spec -sub 40 -num 36 -gain
4 -poly1 - $^{152}\text{Eu}$  |tee recal.log
```

- case 2 : for detectors **without** Trapping.cal files (new detectors undamaged):

```
$ RecalEnergy -spe ../../Out/$CryId/Post__5-40-16384-UI__Ener.spec -sub 0 -num 36 -gain
4 -poly1 -152Eu |tee recal.log
```

This difference is due to the fact that the segment energy corrected from the trapping file is filled in the spectra of the second library of the file (“-sub 40”). But if there is no neutron damage correction (new detectors), these spectra are empty. Spectra of the first library (“-sub 0”) need thus to be taken for the re-calibration.

2) Extract the re-calibration parameters:

```
$ tail -n 36 recal.log |tee recal_nohead.log ; awk -F' ' '{printf "segm %2.2s %6.3f
%.6f \n",$1,$16,$17}' recal_nohead.log |tee RecalEnergy2.cal
```

3) Adapt the *gen_conf.py* as follow:

```
#####
PostPSAFilter=(
"ActualClass          PostPSAFilter",          # name of the used daughter class
"SaveDataDir          $SAVEDIR/$CRYSTAL",        # normally Out/Data(online)
"EnergyGain           4",                      # channels/keV of the calibrated energy spectra
### command lines to be produced only for the specified crystals
{
'00A' : ("TrappingFile Trapping.cal", "RecalEnergy2 RecalEnergy2.cal"),
}
)
```

Check also that the file: RecalEnergy2.cal is written in the “ExtraFiles” list of the *gen_conf.py*:

```
ExtraFiles={
'CRYSTAL' : "CrystalProducerATCA.conf PreprocessingFilterPSA.conf xinv_1325-1340.cal xdir_1325-
1340.cal Trapping.cal RecalEnergy2.cal",
}
```

4) Core recalibration:

- case 1 : for detectors **with** Trapping.cal files **present** (same reasons as for segments):

```
$ RecalEnergy -spe ../../Out/$CryId/Post__5-40-16384-UI__Ener.spec -sub 79 -num 1 -gain
4 -poly1 -152Eu| tee recalCore.log
```

- case 2 : for detectors **without** Trapping.cal files:

```
$ RecalEnergy -spe ../../Out/$CryId/Post__5-40-16384-UI__Ener.spec -sub 39 -num 1 -gain
4 -poly1 -152Eu| tee recalCore.log
```

5) Extract the re-calibration parameters:

```
$ tail -n 1 recalCore.log |tee recalCore_nohead.log ; awk -F' ' '{printf "\"RecalCC
%6.3f %6.3f\" \n",$16,$17}' recalCore_nohead.log
```

The output of the last command should be equivalent to: *"RecalCC -0.115 1.000401"* and needs to be copied in the *gen_conf.py* as follow:

```
#####
PostPSAFilter=(
"ActualClass      PostPSAFilter",      # name of the used daughter class
"SaveDataDir      $SAVEDIR/$CRYSTAL",  # normally Out/Data(online)
"EnergyGain       4",                  # channels/keV of the calibrated energy spectra
### command lines to be produced only for the specified crystals
{
'00A' : ("RecalCC -0.094 1.000651","TrappingFile Trapping.cal", "RecalEnergy2 RecalEnergy2.cal"),
}
)
```

6) Check the results:

To check that the recalibration has been well applied. A new replay is necessary with the modified *gen_conf.py* file. The file "Post__5-40-16384-UI__Ener.spec" can be analyzed with TkT. The different libraries are:

0. segment energies without correction
1. segment energies after trapping corrections
2. fraction of the core energy without correction
3. fraction of the core energy after trapping corrections
- 4. segment energies after trapping corrections and re-calibration.**

Force Segments to Core

The final step at the post PSA level is to decide of the use, or not, of the ForceSegmentsToCore option. If the core energy resolution is good and the counting rates not too important, it is better to use it. This option is defined for the desired detectors in the *gen_conf.py*, as follow:

```
#####
PostPSAFilter=(
"ActualClass      PostPSAFilter",      # name of the used daughter class
"SaveDataDir      $SAVEDIR/$CRYSTAL",  # normally Out/Data(online)
"EnergyGain       4",                  # channels/keV of the calibrated energy spectra
### command lines to be produced only for the specified crystals
{
'00A' : ("RecalCC -0.094 1.000651","TrappingFile Trapping.cal", "RecalEnergy2
RecalEnergy2.cal", "ForceSegmentsToCore"),
}
)
```

Global time alignments

The final time alignment of the cores is done using the output of the tracking. For this, this kind of topology should be used:

```
LOOP CRY 00B 00C 01B 01C 04B 04C 05A 05B 05C 06A 06B 06C

Chain 3      CRY
```

```

Producer      BasicAFP
Filter        PostPSAFilter
Dispatcher    EventBuilder
ENDLOOP

Chain 3       Builder/
Builder       EventBuilder
Filter        TrackingFilterOFT
Consumer      BasicAFC

```

Then, it is necessary to define in the *gen_conf.py* file the mapping of the detectors in the TrackingFilter actor as follow:

```

#####
TrackingFilter=(
"ActualClass      PTrackingFilterOFT",      # name of the used daughter class
"SaveDataDir      $SAVEDIR/$MERGER",        # Out/Merger
"EnergyGain       4",                       # channels/keV of the calibrated energy spectra
"OftParams        0.05 0.02 0.8",           # minprobtrack minprobsing sigma_thet (0==default)
"SourcePosition   0 0 0",                  # source pos with respect to the center of AGATA
"NumGeDets 24",                               # for printing of TkT spectra
"SpecMap 0 0",                                # detID --> specID
"SpecMap 1 1",
"SpecMap 2 2",
"SpecMap 6 3",
...
...
"SpecMap 39 21",
"SpecMap 40 22",
"SpecMap 41 23",
)

```

Note: You need to be in the global folder (where Out, Data, Conf are located)

1) Fit the time spectra:

The tacking actor produces a spectra file named: “Track__N-N-1000-UI__TT.spec”, with N the number of detectors defined in the TrackingFilter part of the *gen_conf.py* file (in the following example, N=24). This file is used for the global time alignment of the cores. The RecalEnergy code is here again used:

```

$ N=24; NN=$(echo "$N*$N" |bc ); RecalEnergy -spe Out/Global/Track__${N}-${N}-1000-
UI__TT.spec -T 500 -num ${NN} |tee recalT.dat ; tail -n ${NN} recalT.dat |tee
recalT_nohead.dat

```

2) Apply the *SolveTT.py* script (this script is located in the agapro/zUseful folder):

```

$ python solveTT.py -f recalT_nohead.dat -n ${N} -c 13 -p 500

```

The end of the output of this script should be similar to:

```

Shifts that minimize Chi2
-0.220
3.565
3.712
1.187
3.751
...
...
3.505

```



```
1.677
-3.947
1.248
-0.931
```

```
Initial:    Average of 506 nonzero values is -499.98900    Chi2 = 692.93223
Corrected:  Average of 506 nonzero values is -499.98900    Chi2 = 17.45832
```

2) Extract the time shift values:

```
$ Nhead=$(echo "$N+3" |bc ); tail -n ${Nhead} solveTT.dat |tee solveTT_tmp.dat |head -n
${N} |tee solveTT.dat ; awk -F' ' '{printf "\"TimeShiftCC %7.3f \"\n",$1}' solveTT.dat
```

The output of this last command give the time shifts to apply in the PostPSA actor for each crystal in the *gen_conf.py* file as follow:

```
#####
PostPSAFilter=(
  "ActualClass      PostPSAFilter",          # name of the used daughter class
  "SaveDataDir      $SAVEDIR/$CRYSTAL",      # normally Out/Data(online)
  "EnergyGain       4",                      # channels/keV of the calibrated energy spectra
  ##### command lines to be produced only for the specified crystals
  {
    '00A' : ("TimeShiftCC -10.859" ,"RecalCC -0.094 1.000651","TrappingFile Trapping.cal",
    "RecalEnergy2 RecalEnergy2.cal","ForceSegmentsToCore"),
  }
)
```

3) After a final global replay, check that the time alignment is correct. For this, open the spectra file “Out/Global/Track__\${N}-\${N}-1000-UI__TT.spec” (in the grid mode). The first set of N spectra will represent the timing between det 0 and the N-1 other... etc. All should be align at 500.